



Electronic Health Information (EHI) Export format is a CSV file.

Single-Patient Electronic Health Information (EHI) Export Instructions for BD HealthSight™ Clinical Advisor and Infection Advisor:

1. Log into HealthSight Clinical Advisor portal (<https://vsi.medmined.com/phoenix>)
2. Create a new patient data report using 'Patient' data type.
3. Run the report, filter the results for specific Medical Record Number (MRN) and download the result using 'Export to Excel'.
4. To get lab results for a patient, create a new report using 'Laboratory and Vitals' data type.
5. Run the report, filter the results for specific MRN and download the results using 'Export to Excel'.
6. To get micro results for a patient, create a new report using 'Microbiology' data type.
7. Run the report, filter the results for specific MRN and download the results using 'Export to Excel'.
8. To get prescription orders for a patient, create a new report using 'Inpatient Medication Orders' data type.
9. Run the report, filter the results for specific MRN and download the results using 'Export to Excel'.
10. To get related surgeries for a patient, create a new report using 'Surgical Procedures' data type.
11. Run the report, filter the results for specific MRN and download the results using 'Export to Excel'.
12. To get device specific information for a patient, create a new report using 'Devices' data type.
13. Run the report, filter the results for specific MRN and download the results using 'Export to Excel'.
14. To save any report in CSV format, open downloaded excel file, click on 'Save a copy' and select 'CSV UTF-8(Comma delimited)(*.csv)' format to save the file in CSV format.

Patient Population Electronic Health Information (EHI) Export Instructions for BD HealthSight™ Clinical Advisor and Infection Advisor:

1. Contact your BD account manager to submit a request to BD Technical Support (TS) to generate this CSV file